

Fast P

Technical Data



Fast P | Features and Accessories

Audiology	
	A 1
Signal processing (channels) / Gain/MPO (handles)	4/—
Hearing programs (pre-configured sound profiles)	3
SpeechMaster	
HD Music (presets)	
TwinPhone ¹⁾	
EchoShield	
Wireless CROS/BICROS ²⁾	
Directionality (channels)	
Narrow Directionality ¹⁾	_
Directional microphone	_
Spatial SpeechFocus ¹⁾	
SpeechFocus	
TruEar™	
Frequency compression	
Extended bandwidth	
Feedback cancellation	•
eWindScreen binaural¹)	
eWindScreen™ (steps)	
Noise Reduction (channels / steps)	_
Speech and noise management (steps)	_
SoundSmoothing [™] (steps)	_
Directional speech enhancement (steps)	_
Adaptive streaming volume ³⁾	_
SoundBrilliance™ 3)	_
Sound equalizer (classes)	_
Spatial Configurator ¹⁾	_
Span ⁴⁾	_
Direction ⁵⁾	
SoundBalance	_
Fitting	
Insitugram	_
Learning (classes) / Data logging	<u> </u>
Acclimatization manager	_
Tinnitus	
Tinnitus noiser	
Static therapy signal (handles / presets)	<u> </u>
Ocean Waves therapy signal (presets)	-

Fast P | Features and Accessories

Style Specific Features	
Ingress Protection Rating	_
Telecoil	_
AutoPhone™	_
Charging contacts	_
Battery Size	13
Battery door on/off function	•
Nanocoated housing	•
e2e wireless™ 3.0	_
Audio streaming	_
User controls coupling via e2e	_
Wireless programming via ConnexxLink™	_
Instrument configurations	
Flat cover	_
Volume wheel	_
Push button	_
Rocker switch	•
Color conversion kit	_
Battery door – direct audio input	_
Battery door – child lock	_

Accessories	
miniPocket	_
CROS Pure	_
eCharger	_
easyPocket™	_
easyTek	
TV Transmitter (req. easyTek)	
Transmitter (req. easyTek)	
VoiceLink™ (req. easyTek)	_
Арр	
MobileFitting App	
easyTek App (req. easyTek)	_
touchControl App	_

lacktriangle available lacktriangle highest feature performance lacktriangle optional lacktriangle not available

 $^{^{\}text{1)}}$ req. bilateral fitting and e2e $^{\text{TM}}$ 3.0

²⁾ req. CROS Pure accessory

 $^{^{\}rm 3)}$ streaming only, req. easyTek $^{\rm TM}$

⁴⁾ req. easyTek & easyTek App, touchControl App or rocker switch

⁵⁾ req. easyTek & easyTek App or touchControl App