



Motion C&G P X

Technical Data

Made for
 iPhone | iPad | iPod

7X

5X

3X

2X

1X

DX



Motion C&G P X | Features and Accessories

| | 7X | 5X | 3X | 2X | 1X |
|---|----------|----------|----------|----------|--------|
| Dynamic Soundscape Processing | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ |
| OVP (Own Voice Processing) ¹⁾ | ■■■■■ | ■■■■■ | ■■■■■ | — | — |
| Sound Clarity | ■■■■■ | ■■■■■ | ■■■ | ■■■ | ■■ |
| Signal processing (channels) / Gain&MPO (handles) | 48 / 20 | 32 / 16 | 24 / 12 | 16 / 8 | 16 / 8 |
| Hearing programs | 6 | 6 | 6 | 4 | 4 |
| Extended dynamic range | ✓ | ✓ | ✓ | ✓ | ✓ |
| Extended bandwidth | ✓ | — | — | — | — |
| EchoShield | ✓ | — | — | — | — |
| HD Music (presets) | 3 | 3 | 1 | 1 | — |
| eWindScreen ²⁾ | Binaural | Binaural | Monaural | Monaural | — |
| Speech and noise management | ✓ | ✓ | ✓ | ✓ | ✓ |
| SoundSmoothing | ✓ | ✓ | ✓ | ✓ | — |
| Feedback cancellation | ✓ | ✓ | ✓ | ✓ | ✓ |
| Speech Quality | ■■■■■ | ■■■■■ | ■■■ | ■■ | ■■ |
| Directionality (Automatic / Adaptive) | Binaural | Binaural | Binaural | ✓ | ✓ |
| Spatial SpeechFocus ^{1) 3)} | ✓ | ✓ | — | — | — |
| TwinPhone ¹⁾ | ✓ | ✓ | ✓ | — | — |
| Frequency compression | ✓ | ✓ | ✓ | ✓ | ✓ |
| Wearer Interaction | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ |
| Signia App (iOS and Android) | ✓ | ✓ | ✓ | ✓ | ✓ |
| Spatial Configurator | ✓ | ✓ | — | — | — |
| Adaptive Streaming Volume ⁴⁾ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Direct Streaming | ✓ | ✓ | ✓ | ✓ | ✓ |
| Made for iPhone | ✓ | ✓ | ✓ | ✓ | ✓ |
| Tinnitus | ✓ | ✓ | ✓ | ✓ | — |
| Notched Amplification Therapy | ✓ | ✓ | ✓ | ✓ | — |
| Tinnitus noise therapy signal | ✓ | ✓ | ✓ | ✓ | — |
| Fitting | ✓ | ✓ | ✓ | ✓ | ✓ |
| Smart Optimizer and Data Logging | ✓ | ✓ | ✓ | ✓ | ✓ |
| Acclimatization manager | ✓ | ✓ | ✓ | ✓ | ✓ |
| InSituGram | ✓ | ✓ | ✓ | ✓ | ✓ |
| AutoFit | ✓ | ✓ | ✓ | ✓ | ✓ |
| TeleCare | ✓ | ✓ | ✓ | ✓ | ✓ |
| Remote Services | ✓ | ✓ | ✓ | ✓ | ✓ |
| Signia App | ✓ | ✓ | ✓ | ✓ | ✓ |

¹⁾ req. bilateral fitting

²⁾ Binaural used in dedicated programs for 5X

³⁾ for 5X, right / left directionality available only in Stroll Program and via the Spatial Configurator

⁴⁾ streaming only

■■■■■ highest feature performance

✓ available — not available ○ optional

Motion C&G P X | Features and Accessories

| | 7X / 5X / 3X | 2X / 1X |
|-----------------------------------|--------------|-----------|
| Style specific features | | |
| Ingress Protection Rating | IP68 | IP68 |
| Charging contacts | ✓ | ✓ |
| Battery Size | — | — |
| Battery door on/off function | — | — |
| Nanocoated housing | ✓ | ✓ |
| e2e wireless 3.0 | ✓ | ✓ |
| User controls coupling via e2e | ✓ | ✓ |
| Wireless programming | ✓ | ✓ |
| Instrument configurations | | |
| Flat cover | — | — |
| Rotary volume control | — | — |
| Push button | — | — |
| Rocker switch | ✓ | ✓ |
| Color conversion kit | ○ | ○ |
| Color conversion kit with T-Coil | — | — |
| T-Coil | ✓ | ✓ |
| Battery door – child lock | — | — |
| Small earhook | ○ | ○ |
| Accessories | | |
| D&C Charger BTE P / Charger BTE P | Mandatory | Mandatory |
| miniPocket | ○ | ○ |
| StreamLine TV | ○ | ○ |
| StreamLine Mic | ○ | ○ |
| CROS Pure 312 X | ○ | — |
| CROS Pure Charge&Go X | ○ | — |
| CROS Silk X | — | — |

✓ available — not available ○ optional



Motion C&G SP X

Technical Data

Made for
 iPhone | iPad | iPod

7X

5X

3X

2X

1X

DX



Motion C&G SP X | Features and Accessories

| | 7X | 5X | 3X | 2X | 1X |
|---|----------|----------|----------|----------|--------|
| Dynamic Soundscape Processing | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ |
| OVP (Own Voice Processing) ¹⁾ | ■■■■■ | ■■■■■ | ■■■■■ | — | — |
| Sound Clarity | ■■■■■ | ■■■■■ | ■■■ | ■■■ | ■■ |
| Signal processing (channels) / Gain&MPO (handles) | 34 / 16 | 32 / 16 | 24 / 12 | 16 / 8 | 16 / 8 |
| Hearing programs | 6 | 6 | 6 | 4 | 4 |
| Extended dynamic range | ✓ | ✓ | ✓ | ✓ | ✓ |
| Extended bandwidth | — | — | — | — | — |
| EchoShield | ✓ | — | — | — | — |
| HD Music (presets) | 3 | 3 | 1 | 1 | — |
| eWindScreen ²⁾ | Binaural | Binaural | Monaural | Monaural | — |
| Speech and noise management | ✓ | ✓ | ✓ | ✓ | ✓ |
| SoundSmoothing | ✓ | ✓ | ✓ | ✓ | — |
| Feedback cancellation | ✓ | ✓ | ✓ | ✓ | ✓ |
| Speech Quality | ■■■■■ | ■■■■■ | ■■■ | ■■ | ■■ |
| Directionality (Automatic / Adaptive) | Binaural | Binaural | Binaural | ✓ | ✓ |
| Spatial SpeechFocus ^{1) 3)} | ✓ | ✓ | — | — | — |
| TwinPhone ¹⁾ | ✓ | ✓ | ✓ | — | — |
| Frequency compression | ✓ | ✓ | ✓ | ✓ | ✓ |
| Wearer Interaction | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ |
| Signia App (iOS and Android) | ✓ | ✓ | ✓ | ✓ | ✓ |
| Spatial Configurator | ✓ | ✓ | — | — | — |
| Adaptive Streaming Volume ⁴⁾ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Direct Streaming | ✓ | ✓ | ✓ | ✓ | ✓ |
| Made for iPhone | ✓ | ✓ | ✓ | ✓ | ✓ |
| Tinnitus | ✓ | ✓ | ✓ | ✓ | — |
| Notched Amplification Therapy | — | — | — | — | — |
| Tinnitus noise therapy signal | ✓ | ✓ | ✓ | ✓ | — |
| Fitting | ✓ | ✓ | ✓ | ✓ | ✓ |
| Smart Optimizer and Data Logging | ✓ | ✓ | ✓ | ✓ | ✓ |
| Acclimatization manager | ✓ | ✓ | ✓ | ✓ | ✓ |
| InSituGram | ✓ | ✓ | ✓ | ✓ | ✓ |
| AutoFit | ✓ | ✓ | ✓ | ✓ | ✓ |
| TeleCare | ✓ | ✓ | ✓ | ✓ | ✓ |
| Remote Services | ✓ | ✓ | ✓ | ✓ | ✓ |
| Signia App | ✓ | ✓ | ✓ | ✓ | ✓ |

¹⁾ req. bilateral fitting

²⁾ Binaural used in dedicated programs for 5X

³⁾ for 5X, right / left directionality available only in Stroll Program and via the Spatial Configurator

⁴⁾ streaming only

■■■■■ highest feature performance

✓ available — not available ○ optional

Motion C&G SP X | Features and Accessories

| | 7X / 5X / 3X | 2X / 1X |
|-------------------------------------|--------------|-----------|
| Style specific features | | |
| Ingress Protection Rating | IP68 | IP68 |
| Charging contacts | ✓ | ✓ |
| Battery Size | — | — |
| Battery door on/off function | — | — |
| Nanocoated housing | ✓ | ✓ |
| e2e wireless 3.0 | ✓ | ✓ |
| User controls coupling via e2e | ✓ | ✓ |
| Wireless programming | ✓ | ✓ |
| Instrument configurations | | |
| Flat cover | — | — |
| Rotary volume control | — | — |
| Push button | — | — |
| Rocker switch | ✓ | ✓ |
| Color conversion kit | ○ | ○ |
| Color conversion kit with T-Coil | — | — |
| T-Coil | ✓ | ✓ |
| Battery door – child lock | — | — |
| Small earhook | ○ | ○ |
| Accessories | | |
| D&C Charger BTE SP / Charger BTE SP | Mandatory | Mandatory |
| miniPocket | ○ | ○ |
| StreamLine TV | ○ | ○ |
| StreamLine Mic | ○ | ○ |
| CROS Pure 312 X | ○ | — |
| CROS Pure Charge&Go X | ○ | — |
| CROS Silk X | — | — |

✓ available — not available ○ optional



Motion C&G X

Technical Data

Made for
 iPhone | iPad | iPod

- 7X
- 5X
- 3X
- 2X
- 1X
- DX



Motion C&G X | Features and Accessories

| | 7X | 5X | 3X | 2X | 1X |
|---|----------|----------|----------|----------|--------|
| Dynamic Soundscape Processing | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ |
| OVP (Own Voice Processing) ¹⁾ | ■■■■■ | ■■■■■ | ■■■■■ | — | — |
| Sound Clarity | ■■■■■ | ■■■■■ | ■■■ | ■■■ | ■■ |
| Signal processing (channels) / Gain&MPO (handles) | 48 / 20 | 32 / 16 | 24 / 12 | 16 / 8 | 16 / 8 |
| Hearing programs | 6 | 6 | 6 | 4 | 4 |
| Extended dynamic range | ✓ | ✓ | ✓ | ✓ | ✓ |
| Extended bandwidth | ✓ | — | — | — | — |
| EchoShield | ✓ | — | — | — | — |
| HD Music (presets) | 3 | 3 | 1 | 1 | — |
| eWindScreen ²⁾ | Binaural | Binaural | Monaural | Monaural | — |
| Speech and noise management | ✓ | ✓ | ✓ | ✓ | ✓ |
| SoundSmoothing | ✓ | ✓ | ✓ | ✓ | — |
| Feedback cancellation | ✓ | ✓ | ✓ | ✓ | ✓ |
| Speech Quality | ■■■■■ | ■■■■■ | ■■■ | ■■ | ■■ |
| Directionality (Automatic / Adaptive) | Binaural | Binaural | Binaural | ✓ | ✓ |
| Spatial SpeechFocus ^{1) 3)} | ✓ | ✓ | — | — | — |
| TwinPhone ¹⁾ | ✓ | ✓ | ✓ | — | — |
| Frequency compression | ✓ | ✓ | ✓ | ✓ | ✓ |
| Wearer Interaction | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ | ■■■■■ |
| Signia App (iOS and Android) | ✓ | ✓ | ✓ | ✓ | ✓ |
| Spatial Configurator | ✓ | ✓ | — | — | — |
| Adaptive Streaming Volume ⁴⁾ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Direct Streaming | ✓ | ✓ | ✓ | ✓ | ✓ |
| Made for iPhone | ✓ | ✓ | ✓ | ✓ | ✓ |
| Tinnitus | ✓ | ✓ | ✓ | ✓ | — |
| Notched Amplification Therapy | ✓ | ✓ | ✓ | ✓ | — |
| Tinnitus noise therapy signal | ✓ | ✓ | ✓ | ✓ | — |
| Fitting | ✓ | ✓ | ✓ | ✓ | ✓ |
| Smart Optimizer and Data Logging | ✓ | ✓ | ✓ | ✓ | ✓ |
| Acclimatization manager | ✓ | ✓ | ✓ | ✓ | ✓ |
| InSituGram | ✓ | ✓ | ✓ | ✓ | ✓ |
| AutoFit | ✓ | ✓ | ✓ | ✓ | ✓ |
| TeleCare | ✓ | ✓ | ✓ | ✓ | ✓ |
| Remote Services | ✓ | ✓ | ✓ | ✓ | ✓ |
| Signia App | ✓ | ✓ | ✓ | ✓ | ✓ |

¹⁾ req. bilateral fitting

²⁾ Binaural used in dedicated programs for 5X

³⁾ for 5X, right / left directionality available only in Stroll Program and via the Spatial Configurator

⁴⁾ streaming only

■■■■■ highest feature performance
 ✓ available — not available ○ optional

Motion C&G X | Features and Accessories

| | 7X / 5X / 3X | 2X / 1X |
|----------------------------------|--------------|-----------|
| Style specific features | | |
| Ingress Protection Rating | IP68 | IP68 |
| Charging contacts | — | — |
| Battery Size | — | — |
| Battery door on/off function | — | — |
| Nanocoated housing | ✓ | ✓ |
| e2e wireless 3.0 | ✓ | ✓ |
| User controls coupling via e2e | ✓ | ✓ |
| Wireless programming | ✓ | ✓ |
| Instrument configurations | | |
| Flat cover | — | — |
| Rotary volume control | — | — |
| Push button | — | — |
| Rocker switch | ✓ | ✓ |
| Color conversion kit | ○ | ○ |
| Color conversion kit with T-Coil | — | — |
| T-Coil | — | — |
| Battery door – child lock | — | — |
| Small earhook | ○ | ○ |
| Accessories | | |
| Inductive Charger II | Mandatory | Mandatory |
| miniPocket | ○ | ○ |
| StreamLine TV | ○ | ○ |
| StreamLine Mic | ○ | ○ |
| CROS Pure 312 X | ○ | — |
| CROS Pure Charge&Go X | ○ | — |
| CROS Silk X | — | — |

✓ available — not available ○ optional